

Fig. 1A






Fig. 1B



Standard curve generated with a cysteine-biotinylated peptide and quantitation by stable isotope dilution. A) Zoom-scan from an ion-trap mass spectrometer showing a 4 amu isotope distribution for the  $[M+2H]^{2+}$  ions of the peptide modified with the isotopically light (1457.9 u) and heavy (1461.8) biotinylating reagents. The ratio (d0/d8) was 4.54. B) Curve generated from the analysis of isotope ratios from zoom-scans of 5 different concentrations of d0-labeled peptide measured in the presence of a known amount of peptide labeled with the isotopically heavy reagent.



Tandem mass spectrum of a cysteine-modified peptide from  $\alpha$ -lactalbumin. Modification of the cysteine residue with the custom synthesized biotinylating reagent did not affect the ability of the Sequest computer program to correctly match this peptide to a database sequence.

Fig. 2

FO00210-4886E850

Fig. 3A




Fig. 3B





Fig. 3C






Fig. 4B


s0319\_hlnavcid.0364.0364.2.out

# amino acids - 93009033, # proteins - 290043, # matched peptides - 1973750

C:\LCQ\database\owl.v31.3, (C# +494.50)


#	Rank/Sp	(M+H) + C*10 <sup>4</sup>	Ions	Reference	Peptide
1.	1 / 1	1994.3 4.4675	17/26	G3P_RABIT	(R)VPTPNVSVVDLTC#R (SEQ ID NO:60)
2.	2 / 403	1995.1 2.7366	13/34	SLTRNGL	(E)LGKPVLTANQVTIWEGLR (SEQ ID NO:61)
3.	3 / 3	1995.0 2.6591	16/36	FLP_LACCA	(N)IANPNVYTETLTAATVCTI (SEQ ID NO:62)
4.	4 / 209	1995.0 2.6335	14/36	A42912	(Y)LALLPSDAEGPHGQFVTDK (SEQ ID NO:63)
5.	5 / 381	1995.1 2.4634	13/38	H69373	(L)ALLVLVAPAMAAGNGEDLRN (SEQ ID NO:64)

Fig. 5A




ADH1 : YSVC#HTDLHAWHGDWPLPVK

ADH2 : YSVC#HTDLHAWHGDWPLPIK



Ratio: 0.57

Fig. 5B



Ratio: >200

Fig. 5C




Figure 6A-C




Figure 7 Schematic representation of the automated LC-MS/MS system. Proteins are typically separated by 1D or 2D SDS-PAGE (1). Protein spots or bands are selected, excised and proteolytically cleaved with trypsin (2). Digests are loaded into an autosampler, which delivers them sequentially to the injection mechanism of a narrow-bore HPLC system (Michrom). The column gradient is automatically applied to separate individual peptides (3). Column eluate is sprayed directly into a mass spectrometer where sequence information from the peptides is collected (4). Recorded peptide masses and CID spectra are transferred to a data station for Sequent analysis, and a final summary of all identifications made for all samples originally loaded is sent to a printer (5).




Figure 8 Schematic of the SPE-CE-MS/MS system. A) A fused silica capillary, typically of 50  $\mu\text{m}$  i.d., is modified at the electrospray end with a liquid junction to establish electrical contact with the analytes inside the capillary. Approx. 5 cm from the end of the capillary, the SPE device is introduced. This consists of  $\text{C}_{18}$ -derivatized, large pore silica beads packed inside a 250  $\mu\text{m}$  i.d. Teflon tubing with Teflon membranes at each end to hold the beads in place between the two fused silica capillaries. B) The injection end of the capillary is inserted into a sealed container which is maintained at a constant, slightly hyperbaric pressure in order to ensure constant flow. A platinum electrode is inserted through the cap, into the container, in order to allow the electrical contact to be made.